

---

# MeaVis Documentation

*Release 0.2.2.dev2*

**Hofheinz' group**

**Dec 28, 2020**



## REFERENCE:

<b>1 Basic example</b>	<b>3</b>
<b>Python Module Index</b>	<b>17</b>
<b>Index</b>	<b>19</b>



MeaVis is a python framework intended to define how **M**easurements have to be run and a programming interface to **V**isualise resulting datasets.

See more details at [ReadTheDocs.io](https://readthedocs.io).



## BASIC EXAMPLE

### 1.1 General configuration of instruments

Let's assume two basic drivers as follow:

```
import vx111

class AgilentB596X(vx111.Instrument):
    def __init__(self, host):
        super().__init__(host=host)

        self.channel = None
        self.mode = None

    def conf():
        if self.channel and self.mode:
            self.write(":SOUR{}:FUNC:MODE {}".format(self.channel, self.mode))

    def output(self, value):
        self.conf()
        if self.channel:
            self.write(
                ":OUTP{} {}".format(self.channel, "ON" if value else "OFF")
            )

    def set_channel(self, channel):
        self.channel = channel
        self.write(":SOUR{}:FUNC:SHAP DC".format(self.channel))
        self.conf()

    def set_mode(self, mode):
        self.mode = mode
        self.conf()

    def set_value(self, value):
        self.write(":SOUR{}:{} {}".format(self.channel, self.mode, value))

class Keysight344XX(vx111.Instrument):
    def __init__(self, host):
        super().__init__(host=host)

        self.write("*CLS")
```

(continues on next page)

(continued from previous page)

```

self.write("*RST")

self.write("TRIG:SOUR IMM")

self.write("CALC:FUNC AVER")
self.write("CALC:STAT ON")

self.ACorDC = None
self.mode = None

def calc_average(self):
    self.write("*CALC:AVER:AVER")
    return float(self.read())

def conf(self):
    if self.ACorDC and self.mode:
        self.write("CONF: {}: {}".format(self.mode, self.ACorDC))

def conf_ACorDC(self, ACorDC):
    self.ACorDC = ACorDC
    self.conf()

def conf_mode(self, mode):
    self.mode = mode
    self.conf()

def count(self, value):
    self.write("SAMP:COUN {}".format(value))

def initiate(self):
    self.write("INIT")

def opc(self):
    self.write("*OPC")

def set_aperture(self, value):
    self.write("{}:APER {}".format(self.mode, value))

```

As you see, this driver is very close to a one-to-one correspondance between the SPCI commands and the methods. Of course, methods can be more complexe, however a simplest driver as possible allow more flexibility.

Then classes specific to MeaVis have to be written to decribe how to use this driver:

- *Drivers*: Front panel instrument.
- *MeaVis classes*: How an experimentalist use the front panel.

```

import drivers

import meavis.tags

@meavis.tags.initialiser("power_source.current_source", mode="CURR")
@meavis.tags.initialiser("power_source.voltage_source", mode="VOLT")
class InitialiserB596X:
    def __init__(self, mode):
        self.mode = mode

```

(continues on next page)

(continued from previous page)

```

def initialise(self, handler, channel):
    handler_channel = drivers.AgilentB596X(**handler)

    handler_channel.set_channel(channel)
    handler_channel.set_mode(self.mode)
    handler_channel.output(True)

    return handler_channel

@meavis.tags.initialiser("multimeter.ac_current_meter")
@meavis.tags.kwargs(mode="CURR", ACorDC="AC")
@meavis.tags.initialiser("multimeter.ac_volt_meter")
@meavis.tags.kwargs(mode="VOLT", ACorDC="AC")
@meavis.tags.initialiser("multimeter.dc_current_meter")
@meavis.tags.kwargs(mode="CURR", ACorDC="DC")
@meavis.tags.initialiser("multimeter.dc_volt_meter")
@meavis.tags.kwargs(mode="VOLT", ACorDC="DC")
class Initialiser344XX:
    def __init__(self, mode, ACorDC):
        self.ACorDC = ACorDC
        self.mode = mode

    def initialise(self, handler, channel):
        handler_channel = drivers.KeySight344XX(**handler)

        handler_channel.conf_mode(self.mode)
        handler_channel.conf_ACorDC(True)

        return handler_channel

@meavis.tags.parameter("power_source.current_source.current")
@meavis.tags.attributes(unit="A", delay=0.1)
@meavis.tags.parameter("power_source.voltage_source.voltage")
@meavis.tags.attributes(unit="V", delay=0.1)
class SourceValue:
    def __init__(self, data):
        self.data = data

    def apply(self, handler, value):
        handler.set_value(value)

@meavis.tags.parameter("multimeter.~.aperture")
@meavis.tags.attributes(unit="s")
class DMMAperture:
    def __init__(self, data):
        self.data = data

    def apply(self, handler, value):
        handler.set_aperture(value)

@meavis.tags.parameter("multimeter.~.average_count")
class DMMCount:
    def __init__(self, data):

```

(continues on next page)

(continued from previous page)

```

        self.data = data

    def apply(self, handler, value):
        handler.count(value)

@meavis.tags.measurement("multimeter.ac_current_meter|dc_current_meter.current")
@meavis.tags.measurement("multimeter.ac_volt_meter|dc_volt_meter.voltage")
class DMMAverage:
    def trigger(self, handler):
        handler.initiate()

    def wait(self, handler):
        handler.opc()
        handler.calc_average()

```

The elements mapped after *kwargs* will be used to initialise the corresponding instrument when required. For example, if a source is used as a voltage source, the statement `handler = meavis_user.InitialiserB596X(mode="VOLT").initialise(/* */) will be executed.`

This file have to be loaded as follow:

```
meavis.instruments.inject(meavis_user._meavis_instruments)
```

Note that it cannot be loaded multiple time, otherwise name collisions will happen.

Up to now the configuration is independant of what we want to measure: it only describes how to use instruments, but not how they are connect or what we want to do.

## 1.2 Experiment-specific configuration of intruments

First we describe how instruments are wired and for which purpose with a YAML file:

```

junction_bias:
  instrument: power_source
  usage: voltage_source
  kwargs:
    addr: 192.168.0.0
  attributes:
    channel: 1
junction_current:
  instrument: multimeter
  usage: dc_current_meter
  kwargs:
    host: 192.168.0.1
  attributes:
    channel: 1
junction_voltage:
  instrument: multimeter
  usage: dc_volt_meter
  kwargs:
    host: 192.168.0.1
  attributes:
    channel: 2

```

The elements mapped after *kwargs* will be used to construct the corresponding instrument when required. For example, the multimeter to measure the junction voltage is constructed with the statement: `handler = meavis_user.ConstructorEthernet(host="192.168.0.1").create()` when required. Moreover the attribute `channel: 2` is used for the initialisation `handler = meavis_user.InitialiserB596X( /* */).initialise( /* */ , channel=2)`.

This file have to be loaded as follow:

```
with open("instances.yaml") as file:
    meavis.instruments.register(yaml.safe_load(file))
```

Note that it cannot be loaded multiple time, otherwise name collisions will happen. After this step, parameters and measurements can be accessed as follow:

```
meavis.parameters.junction_current.aperture([10e-3])
meavis.parameters.junction_current.average_count([100])

meavis.parameters.junction_voltage.aperture([100e-3])
meavis.parameters.junction_voltage.average_count([10])
```

Availed parameters and measurements are displayed in the log output:

```
INFO -- Register power_source constructor [90e97748d6ea9cbb434602eb177a91c685701667]
↪{host: 192.168.0.0}.
INFO -- Register voltage_source initialiser {mode: VOLT} for junction_bias.
INFO -- Register voltage as parameter named junction_bias.voltage.
INFO -- Register multimeter constructor [1ca226d3ca09e4167fbbfa3bd218a8323d76e12f]
↪{host: 192.168.0.1}.
INFO -- Register dc_current_meter initialiser {mode: CURR, ACorDC: DC} for junction_
↪current.
INFO -- Register aperture as parameter named junction_current.aperture.
INFO -- Register average_count as parameter named junction_current.average_count.
INFO -- Register current as measurement named junction_current.current.
INFO -- Register dc_volt_meter initialiser {mode: VOLT, ACorDC: DC} for junction_
↪voltage.
INFO -- Register aperture as parameter named junction_voltage.aperture.
INFO -- Register average_count as parameter named junction_voltage.average_count.
INFO -- Register voltage as measurement named junction_voltage.voltage.
INFO -- Add completer group for junction_bias : {junction_bias.voltage}.
INFO -- Add completer group for junction_current : {junction_current.average_count, ↪
↪junction_current.aperture}.
INFO -- Add completer group for junction_voltage : {junction_voltage.aperture, ↪
↪junction_voltage.average_count}.
```

And finally the measurement can be described and processed as follow:

```
measurement_loop = meavis.loop.LoopEngine(
    yaml.safe_load(
        """
parameters:
  - junction_bias.voltage
measurements:
  - junction_current.current
  - junction_voltage.voltage
name: iv_dc_4probes
"""
    )).create(
    meavis.parameters.junction_bias.voltage(numpy.linspace(-1e-3, 1e-3, 401)),
```

(continues on next page)

(continued from previous page)

```

    meavis.measurements.junction_current.current(),
    meavis.measurements.junction_voltage.voltage(),
)
measurement_loop.trigger(None)
measurement_loop.wait(None)

```

In the log output, instruments are created and initialised when required:

```

INFO -- Complete parameters with [junction_current.average_count, junction_current.
↪aperture]
INFO -- Complete parameters with [junction_voltage.aperture, junction_voltage.
↪average_count]
INFO -- Create handler of multimeter.constructor_
↪[1ca226d3ca09e4167fbbfa3bd218a8323d76e12f] with {host: 192.168.0.1}.
INFO -- Initialise channel 2 on handler of multimeter.dc_volt_meter.initialiser with
↪{mode: VOLT, ACorDC: DC}.
INFO -- Set junction_voltage.aperture to 0.1 s.
INFO -- Set junction_voltage.average_count to 10.
INFO -- Initialise channel 1 on handler of multimeter.dc_current_meter.initialiser_
↪with {mode: CURR, ACorDC: DC}.
INFO -- Set junction_current.aperture to 0.01 s.
INFO -- Set junction_current.average_count to 100.
INFO -- Create handler of power_source.constructor_
↪[90e97748d6ea9cbb434602eb177a91c685701667] with {host: 192.168.0.0}.
INFO -- Initialise channel 1 on handler of power_source.voltage_source.initialiser_
↪with {mode: VOLT}.
INFO -- Set junction_bias.voltage to -0.001 V.
INFO -- Trigger junction_current.current, waiting for data.
INFO -- Trigger junction_voltage.voltage, waiting for data.
INFO -- Set junction_bias.voltage to -0.000998 V.
INFO -- Trigger junction_current.current, waiting for data.
INFO -- Trigger junction_voltage.voltage, waiting for data.

```

## 1.2.1 API

### meavis package

Measurement & Visualisation python framework.

### Submodules

#### meavis.completer module

MeaVis parameters completion.

```

class meavis.completer.CompleterEngine(data)
    Bases: object

    Define how parameters of a loop has to be completed.

    instances_parameters = {}

    __init__(data)
        Store a data structure as loop pattern.

```

```
classmethod clear()
    Clear instances set.

classmethod inject_instances(instances)
    Inject instances in CompleterEngine maps.

complete(instances=())
    Complete a loop pattern.
```

## meavis.instruments module

Main loop functions for running MeaVis measurements.

```
meavis.instruments.clear(module_name)
    Clear injected names by users.

meavis.instruments.inject(instruments)
    Inject instruments.

meavis.instruments.register(instances)
    Inject instances.
```

## meavis.loop module

Main loop functions for running MeaVis measurements.

```
class meavis.loop.LoopMeasurement(parameters, measurements)
    Bases: object

    Define a measurement running a loop.

    __init__(parameters, measurements)
        Ceate a loop measurement.

    trigger(handler)
        Nothing to trigger.

    wait(handler)
        Run the loop.

class meavis.loop.LoopEngine(data)
    Bases: object

    Define how a loop has to be processed.

    items_map = {}

    default_map = {'measurements': {'handler': None, 'initialiser': None, 'invasive': None}}

    __init__(data)
        Store a data structure as loop pattern.

    classmethod clear()
        Clear MeaVis item maps.

    classmethod inject_items(*items)
        Inject MeaVis items in the LoopEngine.

    inject_defaults()
        Inject default attributes in LoopEngine maps.
```

**complete()**  
Complete current data structure.

**create** (\*items, completion=True)  
Create a measurement from the pattern.

**synchronisers** (state\_parameters, completion=True)  
Synchronise parameters group from the pattern.

### meavis.markup module

Read and write MeaVis markup language.

**meavis.markup.visit\_instruments** (meavis\_instruments, tag\_name, meavis\_name)  
Visit an instrument hierarchy and return corresponding items.

### meavis.measurements module

MeaVis measurements namespace for user-defined injection.

**meavis.measurements.inject** (cls, name)  
Wrap and inject user-defined measurement in this namespace.

### meavis.parameters module

MeaVis parameters namespace for user-defined injection.

**meavis.parameters.inject** (cls, name)  
Wrap and inject user-defined parameter in this namespace.

### meavis.synchroniser module

Synchronisation of MeaVis parameters.

**class** meavis.synchroniser.**LoopSynchroniser** (state\_parameters, loop\_parameters, synchronisers)  
Bases: `object`  
Define a loop to pre-synchronise instruments.

**\_\_init\_\_** (state\_parameters, loop\_parameters, synchronisers)  
Ceate a loop synchroniser.

**pre\_synchronise** (states)  
Run the loop.

## meavis.tags module

Class decorators to tag MeaVis classes.

`meavis.tags.add_metadata(metadata_name, **kwargs)`

Add metadata to a tagged MeaVis class.

`meavis.tags.add_item(tag_name, meavis_name, **kwargs)`

Tag a MeaVis class.

`meavis.tags.attributes(**kwargs)`

Add attributes to a tagged MeaVis class.

`meavis.tags.kwargs(**kwargs)`

Add kwargs to a tagged MeaVis class.

`meavis.tags.constructor(meavis_name, **kwargs)`

Tag a class as a MeaVis initialiser.

`meavis.tags.initialiser(meavis_name, **kwargs)`

Tag a class as a MeaVis initialiser.

`meavis.tags.measurement(meavis_name, **kwargs)`

Tag a class as a MeaVis measurement.

`meavis.tags.parameter(meavis_name, **kwargs)`

Tag a class as a MeaVis parameter.

## meavis.tasks module

MeaVis tasks for threading.

`meavis.tasks.setup_and_acquire(meavis_item)`

Set meavis\_item handler if required.

`meavis.tasks.settle(parameter, sample, lock_in, lock_out, delay=True)`

Settle a parameter.

`meavis.tasks.trigger_wait(measurement, lock_in, lock_out, lock_barrier)`

Trigger & wait for a measurement.

## 1.2.2 ChangeLog

All notable changes to this project will be documented in this file.

The format is based on [Keep a Changelog](#) and this project adheres to [Semantic Versioning](#).

### Unreleased

#### Added

#### Change

#### Deprecated

## Removed

## Fixed

## [0.2.1]

## Added

- Basic IV measurement as first example.
- Loop completer for parameters.
- Tags mechanism.

## Change

- Mapping inheritance for MeaVis Markup Language – Instrument.
- Constructor and Initialiser tags for MeaVis Markup Language – Instrument.

## 1.2.3 License

GNU LESSER GENERAL PUBLIC LICENSE  
Version 3, 29 June 2007

Copyright (C) 2007 Free Software Foundation, Inc. <<https://fsf.org/>>  
Everyone **is** permitted to copy **and** distribute verbatim copies  
of this license document, but changing it **is not** allowed.

This version of the GNU Lesser General Public License incorporates  
the terms **and** conditions of version 3 of the GNU General Public  
License, supplemented by the additional permissions listed below.

### 0. Additional Definitions.

As used herein, "**this License**" refers to version 3 of the GNU Lesser  
General Public License, **and** the "**GNU GPL**" refers to version 3 of the GNU  
General Public License.

"**The Library**" refers to a covered work governed by this License,  
other than an Application **or** a Combined Work **as** defined below.

An "**Application**" **is** any work that makes use of an interface provided  
by the Library, but which **is not** otherwise based on the Library.  
Defining a subclass of a **class defined** by the Library **is** deemed a mode  
of using an interface provided by the Library.

A "**Combined Work**" **is** a work produced by combining **or** linking an  
Application **with** the Library. The particular version of the Library  
**with** which the Combined Work was made **is** also called the "**Linked  
Version**".

The "**Minimal Corresponding Source**" **for** a Combined Work means the

(continues on next page)

(continued from previous page)

Corresponding Source **for** the Combined Work, excluding **any** source code **for** portions of the Combined Work that, considered **in** isolation, are based on the Application, **and not** on the Linked Version.

The "Corresponding Application Code" **for** a Combined Work means the **object** code **and/or** source code **for** the Application, including **any** data **and** utility programs needed **for** reproducing the Combined Work **from the** Application, but excluding the System Libraries of the Combined Work.

#### 1. Exception to Section 3 of the GNU GPL.

You may convey a covered work under sections 3 **and** 4 of this License without being bound by section 3 of the GNU GPL.

#### 2. Conveying Modified Versions.

If you modify a copy of the Library, **and, in** your modifications, a facility refers to a function **or** data to be supplied by an Application that uses the facility (other than **as** an argument passed when the facility **is** invoked), then you may convey a copy of the modified version:

- a) under this License, provided that you make a good faith effort to ensure that, **in** the event an Application does **not** supply the function **or** data, the facility still operates, **and** performs whatever part of its purpose remains meaningful, **or**
- b) under the GNU GPL, **with** none of the additional permissions of this License applicable to that copy.

#### 3. Object Code Incorporating Material **from Library** Header Files.

The **object** code form of an Application may incorporate material **from** a header file that **is** part of the Library. You may convey such **object** code under terms of your choice, provided that, **if** the incorporated material **is not** limited to numerical parameters, data structure layouts **and** accessors, **or** small macros, inline functions **and** templates (ten **or** fewer lines **in** length), you do both of the following:

- a) Give prominent notice **with** each copy of the **object** code that the Library **is** used **in** it **and** that the Library **and** its use are covered by this License.
- b) Accompany the **object** code **with** a copy of the GNU GPL **and** this license document.

#### 4. Combined Works.

You may convey a Combined Work under terms of your choice that, taken together, effectively do **not** restrict modification of the portions of the Library contained **in** the Combined Work **and** reverse engineering **for** debugging such modifications, **if** you also do each of the following:

- a) Give prominent notice **with** each copy of the Combined Work that the Library **is** used **in** it **and** that the Library **and** its use are covered by this License.

(continues on next page)

(continued from previous page)

b) Accompany the Combined Work **with** a copy of the GNU GPL **and** this license document.

c) For a Combined Work that displays copyright notices during execution, include the copyright notice **for** the Library among these notices, **as** well **as** a reference directing the user to the copies of the GNU GPL **and** this license document.

d) Do one of the following:

0) Convey the Minimal Corresponding Source under the terms of this License, **and** the Corresponding Application Code **in** a form suitable **for**, **and** under terms that permit, the user to recombine **or** relink the Application **with** a modified version of the Linked Version to produce a modified Combined Work, **in** the manner specified by section 6 of the GNU GPL **for** conveying Corresponding Source.

1) Use a suitable shared library mechanism **for** linking **with** the Library. A suitable mechanism **is** one that (a) uses at run time a copy of the Library already present on the user's **computer** system, **and** (b) will operate properly **with** a modified version of the Library that **is** interface-compatible **with** the Linked Version.

e) Provide Installation Information, but only **if** you would otherwise be required to provide such information under section 6 of the GNU GPL, **and** only to the extent that such information **is** necessary to install **and** execute a modified version of the Combined Work produced by recombining **or** relinking the Application **with** a modified version of the Linked Version. (If you use option 4d0, the Installation Information must accompany the Minimal Corresponding Source **and** Corresponding Application Code. If you use option 4d1, you must provide the Installation Information **in** the manner specified by section 6 of the GNU GPL **for** conveying Corresponding Source.)

## 5. Combined Libraries.

You may place library facilities that are a work based on the Library side by side **in** a single library together **with** other library facilities that are **not** Applications **and** are **not** covered by this License, **and** convey such a combined library under terms of your choice, **if** you do both of the following:

a) Accompany the combined library **with** a copy of the same work based on the Library, uncombined **with** any other library facilities, conveyed under the terms of this License.

b) Give prominent notice **with** the combined library that part of it **is** a work based on the Library, **and** explaining where to find the accompanying uncombined form of the same work.

## 6. Revised Versions of the GNU Lesser General Public License.

The Free Software Foundation may publish revised **and/or** new versions

(continues on next page)

(continued from previous page)

of the GNU Lesser General Public License **from time** to time. Such new versions will be similar **in** spirit to the present version, but may differ **in** detail to address new problems **or** concerns.

Each version **is** given a distinguishing version number. If the Library **as** you received it specifies that a certain numbered version of the GNU Lesser General Public License "**or any later version**" applies to it, you have the option of following the terms **and** conditions either of that published version **or** of **any** later version published by the Free Software Foundation. If the Library **as** you received it does **not** specify a version number of the GNU Lesser General Public License, you may choose **any** version of the GNU Lesser General Public License ever published by the Free Software Foundation.

If the Library **as** you received it specifies that a proxy can decide whether future versions of the GNU Lesser General Public License shall apply, that proxy's **public statement of acceptance of any version is** permanent authorization **for** you to choose that version **for** the Library.



## PYTHON MODULE INDEX

### m

- meavis, 8
- meavis.completer, 8
- meavis.instruments, 9
- meavis.loop, 9
- meavis.markup, 10
- meavis.measurements, 10
- meavis.parameters, 10
- meavis.synchroniser, 10
- meavis.tags, 11
- meavis.tasks, 11



## Symbols

`__init__()` (*meavis.completer.CompleterEngine method*), 8  
`__init__()` (*meavis.loop.LoopEngine method*), 9  
`__init__()` (*meavis.loop.LoopMeasurement method*), 9  
`__init__()` (*meavis.synchroniser.LoopSynchroniser method*), 10

## A

`add_item()` (*in module meavis.tags*), 11  
`add_metadata()` (*in module meavis.tags*), 11  
`attributes()` (*in module meavis.tags*), 11

## C

`clear()` (*in module meavis.instruments*), 9  
`clear()` (*meavis.completer.CompleterEngine class method*), 8  
`clear()` (*meavis.loop.LoopEngine class method*), 9  
`complete()` (*meavis.completer.CompleterEngine method*), 9  
`complete()` (*meavis.loop.LoopEngine method*), 9  
`CompleterEngine` (*class in meavis.completer*), 8  
`constructor()` (*in module meavis.tags*), 11  
`create()` (*meavis.loop.LoopEngine method*), 10

## D

`default_map` (*meavis.loop.LoopEngine attribute*), 9

## I

`initialiser()` (*in module meavis.tags*), 11  
`inject()` (*in module meavis.instruments*), 9  
`inject()` (*in module meavis.measurements*), 10  
`inject()` (*in module meavis.parameters*), 10  
`inject_defaults()` (*meavis.loop.LoopEngine method*), 9  
`inject_instances()` (*meavis.completer.CompleterEngine class method*), 9  
`inject_items()` (*meavis.loop.LoopEngine class method*), 9

`instances_parameters` (*meavis.completer.CompleterEngine attribute*), 8  
`items_map` (*meavis.loop.LoopEngine attribute*), 9

## K

`kwargs()` (*in module meavis.tags*), 11

## L

`LoopEngine` (*class in meavis.loop*), 9  
`LoopMeasurement` (*class in meavis.loop*), 9  
`LoopSynchroniser` (*class in meavis.synchroniser*), 10

## M

`measurement()` (*in module meavis.tags*), 11  
`meavis` (*module*), 8  
`meavis.completer` (*module*), 8  
`meavis.instruments` (*module*), 9  
`meavis.loop` (*module*), 9  
`meavis.markup` (*module*), 10  
`meavis.measurements` (*module*), 10  
`meavis.parameters` (*module*), 10  
`meavis.synchroniser` (*module*), 10  
`meavis.tags` (*module*), 11  
`meavis.tasks` (*module*), 11

## P

`parameter()` (*in module meavis.tags*), 11  
`pre_synchronise()` (*meavis.synchroniser.LoopSynchroniser method*), 10

## R

`register()` (*in module meavis.instruments*), 9

## S

`settle()` (*in module meavis.tasks*), 11  
`setup_and_acquire()` (*in module meavis.tasks*), 11  
`synchronisers()` (*meavis.loop.LoopEngine method*), 10

## T

`trigger()` (*meavis.loop.LoopMeasurement method*), [9](#)  
`trigger_wait()` (*in module meavis.tasks*), [11](#)

## V

`visit_instruments()` (*in module meavis.markup*),  
[10](#)

## W

`wait()` (*meavis.loop.LoopMeasurement method*), [9](#)